

***Azadinium* isolation from sediment samples**

Joo-Hwan Kim, Brian Bill, Nicolaus Adams and Vera
L. Trainer

Hanyang University, Seoul, South Korea
NOAA-Fisheries Science Center, Seattle, WA, USA

❖ Background

- ◆ Sediment incubation method has been used to isolate *Azadinium*.
- ◆ According to Adams, cell concentration of *Azadinium* spp. in Puget Sound is low (< 200 cells/L).
- ◆ Low cyst abundance of *Azadinium* spp. is likewise expected in Puget Sound area.


❖ Background

- ◆ It could be time, labor, money consuming work that incubation and observation of all sediment samples.
- ◆ qPCR assay on sediment samples could be helpful method for pre-screening of *Azadinium* cysts.
- ◆ We tried to screening of *Azadinium* cysts using qPCR and selected sediment samples for incubation & cell isolation

❖ Materials & Methods

◆ Sediment treatment

- 1-2 g of sediment was diluted 10 mL of filtered seawater (0.45 μm) to make sediment suspension.




DNA extraction



Sediment incubation

◆ DNA extraction

- We used heating method to remove DNA debris in sediment samples which could seriously overestimation effect on result of qPCR assay (Kim et al. 2016).



◆ qPCR assay

- We used Amphidomataceae and *A. poporum*, *A. obesum*, *A. spinosum* specific primer for qPCR assay (Toebe, 2013; Smith, 2015).

Table II: FISH and qPCR primer and probes specific for Azadinium species

Target species and target gene	Target nucleotide accession numbers	FISH probe (5'–3')	TaqMan MGB probe (5'–3')	Primer forward (5'–3')	Primer reverse (5'–3')	qPCR amplicon size [number of base pairs (bp)]
<i>A. spinosum</i> , 28S	HQ324896	Asp_544	Aspin77T	Asp48F	Asp120R	72 bp
	JN165101	TGG TCG AGT	CGC CCA AAA	TCG TCT TTG TGT CAG	GGA AAC TCC TGA	
	FJ217815	TAC CAG CCC	GGA CTC CT	GGA GAT G	AGG GCT TGT	
<i>A. poporum</i> , 28S	HQ324893	Apop_544	Apop112	Apop62F	Apop148R	68 bp
	HQ324894	CGA GTT ACC	TTC CAG ACG	GAT GCT CAA GGT GCC	CCT GCG TGT CTG	
	HQ324895	AGT TCT CCG	ACT CAA A	TAG AAA GTC	GTT GCA	
<i>A. obesum</i> , 28S	GQ914936	Aob_544	Aob163	Aob134F	Aob208R	74 bp
		AAG ACA TTC	AAG ACA TTC	AGG GAT CGA TAC ACA	AAA CTC CAG GGA	
		GAC CTA CCG	GAC CTA CCG T	AAT GAG TAC TG	CAT GGT AGT CTT A	


Table 1 Sequences of primers and probes used in this study including optimised final concentrations and annealing temperatures for real-time PCR assays

Target	Type	Sequence	Product size	Annealing temperature	Final concentration	Reference
Amphidomataceae						This study
Amp240F	Forward primer	CAA CTT TCA GCG ACG GAT GTC TCG	179 bp	62°C	200 nM	
Amp418R	Reverse primer	AAG CYR CWGGCA TKA GAA CGT ACW CCC			200 nM	



Results

◆ Amphidomataceae qPCR assay on sediment samples (1/16-18/2016)



- Sediment DNA was extracted including DNA debris removal step.


- qPCR assay was conducted using Amphidomataceae specific primer (SYBR-Green method).



Relative abundance of Amphidomataceae cysts

* Cyst abundance value was very roughly calculated based on C_t value of qPCR assay

◆ *Azadinium* specific qPCR assay on sediment samples (1/16-18/2016)





- qPCR assay was conducted using each specific probes of *A. spinosum*, *A. obesum*, *A. poporum*. (Taqman probe method)
- *A. spinosum*, *A. obesum* was not detected in all stations.



Relative abundance of *A. poporum* cysts

* Cyst value was roughly calculated based on C_t value of qPCR assay




n Amphidomataceae cysts

○ *A. poporum* cysts

- The phenomenon of discordance was observed between cyst distributions of family level and species level.
- This data suggests that other *Azadinium* species could exist in Puget Sound area.

◆ Sediment incubation & Isolation

- After qPCR screening, we selected sediment samples for incubation.
- 1-5 mL of sediment suspension was diluted with 5-50 mL of ESNW-Si medium.
- Incubation in 18°C, 14:10 (L:D)
- After 5 days, *Azadinium* like cells were isolated.



Possible *Azadinium* (96-2 b3)

◆ *Azadinium* strains list


- *Azadinium* has been never successfully isolated in Pacific coastal of North America due to their extremely low abundance (detected only PCR results).
- Eight possible *Azadinium* strains were isolated.

Species	Strains	Collection date	DNA extraction	RT PCR assay *
Possible <i>Azadinium</i>	48-1 b2	2/1/2016	-	-
Possible <i>Azadinium</i>	48-1 b3	2/1/2016	-	-
Possible <i>Azadinium</i>	48-1 b5	2/1/2016	-	-
Possible <i>Azadinium</i>	48-1 f2	2/2/2016	o	o
Possible <i>Azadinium</i>	48-1 f8	2/2/2016	o	o
Possible <i>Azadinium</i>	96-2 b3	2/2/2016	o	o
Possible <i>Azadinium</i>	96-2 b8	2/2/2016	o	o
Possible <i>Azadinium</i>	96-1 f4	1/15/2016	-	-


* RT-PCR was conducted using Amphidomataceae specific primer.


◆ Morphological Characteristics



48-1 f2 (*A. obesum*?)



96-2 b3 *A. trinitatum* ?
A. dexteroporum ?







Thank you for your attention !